

Project Name: RAMS2

# Test Case Template

Test Case ID: RAMS2\_T

Test Priority (Low/Medium/High): Med

Module Name: ECGAnalyzer

Test Title: Verify functionality of code for ECGAnalyzer

Description: Testing input, output, and basic functionality, error reporting, and validating basic algorithm

Test Designed by: EmSt

Test Designed date: 07.01.2026

Test Executed by: EmSt

Test Execution date: 07.01.2026

Pre-conditions: User has installed Python Environment and Dependencies

Dependencies: requirements.txt

Step	Test Steps	Test Data	Expected Result	Actual Result	Notes	Status
UT_LC_01	Load valid config and verify attributes + sampling rate	conf.ini	Object created; required attrs exist; samplingRate_hz==200.0	Analyzer initialized; attrs present; samplingRate_hz=200.0	—	Pass
UT_LC_02	Load missing config	/nonexistent/path/conf.ini	FileNotFoundError	FileNotFoundError raised	—	Pass
UT_LC_03	Load invalid config value	temp.ini with samplingRate_hz = not_a_number	ValueError	ValueError raised	Log: “Fehler beim Parsen der Config ... not_a_number”	Pass
UT_LMS_01	Load valid ECG signal from .mat	samples/ecg1.mat	numpy.ndarray, non-empty	Array loaded, non-empty	—	Pass
UT_LMS_02	Load missing .mat file	/nonexistent/file.mat	None returned (or graceful handling)	None returned	Log: “Datei nicht gefunden ...”	Pass
UT_COS_01	Clip outliers in signal	[0.1, 0.2, 0.3, 0.4, 0.5, 10, -10, 0.3, 0.4]	Extremes are clipped (max<10, min>-10)	Outliers clipped	—	Pass
UT_COS_02	Clip outliers on empty array	[]	empty array returned	empty array returned	—	Pass
UT_DPS_01	Detect peaks on synthetic ECG-like signal	length=2000, peaks at [250, 270, 400, ..., 1890]	8 peaks detected (tolerant range)	Peak count in expected range	—	Pass
UT_DPS_02	Detect peaks on empty signal	[]	empty peaks result	empty peaks result	—	Pass
UT_DPS_03	Detect peaks on constant signal	np.ones(1000)*0.5	“No/very few peaks” ( $\leq 2$ )	<b>20 peaks detected</b> at [0, 51, 102, ..., 969]	The peak detector is treating a flat signal as peak-worthy (likely threshold/gradient logic). Either fix algorithm or relax test to check “no physiological peaks”	Fail

					differently.	
UT_CHM_01	HRV metrics for known NN intervals	[ 800, 850, 900, 850, 800 ] ms, duration=60s	HR_Mean≈60000/mean (nn) , SDNN≈std (ddof=1)	Values match within tolerance	—	Pass
UT_CHM_02	HRV metrics with empty NN list	[ ], duration=60s	None returned	None returned	—	Pass
UT_CHM_03	NN50 count correctness	[ 800, 900, 850, 950, 800 ]	NN50=3 (100, 50, 100, 150 → 3 > 50)	NN50=3	—	Pass
UT_CHM_04	SDANN not available for <5min	[ 800 ] *100, duration=60s	SDANN is NaN/N/A	SDANN is NaN	—	Pass
UT_CHM_05	HR min/max correctness	[ 500, 1000, 750 ] ms	HR_Min=60, HR_Max=120	HR_Min/Max match tolerance	—	Pass
ST_AF_01	End-to-end analysis with valid intervals; check printed metrics	intervals [ ( 0, 350 ) , ( 350, 700 ) ]	Output contains Intervall 1/2, HR, SDNN, NN50	Output contains all expected markers	—	Pass
ST_AF_02	>7 intervals triggers summary	8 intervals ( 0..2800 )	Output contains ZUSAMMENFASSUNG, HR Mean, SDNN Mean, Intervall 8	Summary printed; Intervall 8 present	Log warning: “Unphysiologische HR: 25-235 bpm”	Pass
ST_AF_03	Exactly 7 intervals → no summary	7 intervals ( 0..2450 )	No ZUSAMMENFASSUNG; contains Intervall 7	No summary; Intervall 7 present	—	Pass
ST_AF_04	Invalid intervals skipped; valid interval still processed	[ ( 500, 400 ) , ( -10, 100 ) , ( 0, 9999999 ) , ( 100, 200 ) ]	Program doesn’t crash; valid interval processed	Valid interval produced HR output	Logs: invalid intervals aborted (3 warnings)	Pass
ST_AF_05	Very short interval doesn’t crash	[ ( 0, 1 ) ]	No crash (may warn “too few peaks”)	No crash; warning emitted	Log: “Zu wenige Peaks ... Abbruch.”	Pass
ST_AF_06	Long interval allows SDANN	[ ( 0, 400 ) ]	Output contains SDANN and not N/A	SDANN present and not N/A	—	Pass
ST_AF_07	Missing file handled gracefully	file=/nonexistent /file.mat, intervals [ ( 0, 100 ) ]	No crash; graceful error handling	No crash; errors logged	Logs: file not found + “Kein Signal erhalten ...”	Pass
ST_AF_08	Empty interval list handled	intervals [ ]	No crash	No crash	—	Pass
ST_AF_09	HR plausibility check	intervals [ ( 0, 60 ) ]	Extracted HR within 30–220 bpm	HR found and within range	—	Pass
ST_AF_10	Repeatability across runs	intervals [ ( 0, 100 ) ] run twice	HR identical within tolerance	HR1 ≈ HR2 (abs<0.01)	—	Pass
ST_INT_01	Full workflow w/ main.py intervals	8 intervals [ ( 0, 350 ) ... ( 4400, 5300 ) ]	Intervall 1..8 + ZUSAMMENFASSUNG + BERICHT	All markers found	Logs: warnings about unphysiological HR ranges	Pass
ST_EDGE_01 (intervals0)	Edge-case intervals do not crash	[ ( 0.0, 60.0 ) ]	No crash	No crash	Float seconds accepted	Pass
ST_EDGE_01 (intervals1)	Edge-case intervals do not crash	[ ( 0, 60 ) , ( 10, 70 ) ]	No crash	No crash	Overlapping intervals tolerated	Pass
ST_EDGE_01 (intervals2)	Edge-case intervals do not crash	[ ( 0, 0 ) ]	No crash (likely interval rejected)	No crash; interval aborted	Log: “Ungültiges Intervall [0, 0]. Abbruch.”	Pass

**Post-conditions:**

Test results have successfully passed. Launching main.py works with the delivered code.